



UNIVERSITÀ
DEGLI STUDI
DI PADOVA



Lymphoma Subtype Classification Using Neural Networks To Support Human Hematopathologist decisions

Chiara Fantinato, Riccardo Malgarotto Boschiero

Human Data Analytics
Academic Year 2022/2023

- ☐ **Introduction**
- ☐ **Dataset**
- ☐ **Pre-processing**
- ☐ **Processing Pipeline**
 - ☐ **ResNet50-based architecture**
 - ☐ **ResNet built from scratch**
 - ☐ **Heatmaps**
 - ☐ **Visual attention mechanism**
- ☐ **Results**
- ☐ **Concluding Remarks**

Background


- **Lymphoma** is a type of cancer that begins in lymphocytes.
- The delivered **treatment** depends on the type of the lymphoma.
- **Lymphoma subtype classification** is a **complex task** even for expert hematopathologists.
- Implementation of algorithms for the **automated classification** of Lymphoma subtypes could be helpful to support physicians' decisions.
- Deep learning models struggles to show impact in medical domain due to a **lack of transparency**.

Aim of the study


- Perform **lymphoma subtype classification** using neural networks.
- Provide **visual explanations** for machine's decisions.

- 374 images of size 1388x1040
 - 113 images for Chronic Lymphocytic Leukemia (CLL)
 - 139 images for Follicular Lymphoma (FL)
 - 122 images Mantle Cell Lymphoma (MCL)
- Large degree of staining and sectioning **variation**.
- Randomly divided into **training, validation, and test set** according to **6:2:2**.


CLL



FL




MCL




Artifact management

- 3x3 gaussian filter.
- Images affected by artifacts **detected automatically**.
 - **Variance** of image filtered using a **Laplacian filter** greater than 700 → Presence of artifacts


No artifacts




Artifacts



Filtered




Pre-processing (2)




- **Histogram equalization** → improve image contrast
- **Macenko normalization** → reduce the staining variation in histological slides


Original




Histogram eq.



Macenko norm.




Hist eq. + M. norm.



Data Augmentation

- Split images into **patches**
 - Resize images from 1388x1040 to 1300x1040 and split into 20 non-overlapping patches of 260x260 with a stride of 260
 - Resize images from 1388x1040 to 1300x1040 and split into 63 overlapping patches of 260x260 with a stride of 130
 - Probability of belonging to a class is obtained by summing the output probabilities of each patch
- **Rotate and/or flip** images (after **resize** from 1388x1040 to 224x224): 8 rotated and flipped versions of each image were retrieved





Basic elements of ResNets:



- Convolutional block
- Residual block

ResNet50


Transfer learning: weights initialized weights computed on ImageNet and then adapted to our dataset.




Proposed ResNet Architecture





- **Idea:** Provide visual explanation for machine's decisions.
- **Heatmaps:** Weighted average of the last convolution layers using as weights the parameters the last dense layer.
- **Validity of heatmaps:**
 - Mean of the ratios between
 - Number of non-zero pixels in the intersection image
 - Number of non-zero pixels in the image obtained by summing all the masked images related to that image
 - Number of images without intersection




Visual attention mechanism



Idea: Histological lymphoma features could be focused on a restricted area of the slide.



Results (1)




- **ResNet50:** 10 epochs, 32 batch size, AdaMax optimizer, variable learning rate from 0.001 to 0.0001
- **Proposed ResNet:** 15 epochs, 32 batch size, AdaMax optimizer, variable learning rate from 0.001 to 0.0001
- **Best pre-processing:** gaussian filter to remove artifacts, non-overlapping patches data augmentation (histogram equalization and Macenko normalization are not useful)

Model	Test accuracy [%]
Orlov et al. [2]	98.00
Andrew Janowczyk and Anant Madabhushi [4]	96.58
Tambe et al. [5]	97.33
Zhang et al. [6]	98.63
Hathem et al. [7]	98.70
ResNet50	98.67
Proposed ResNet	94.67

Results (2)




- Our ResNet misclassifies 4 images out of 75 in the test set.
- $\text{TPR}_{\text{MCL}} = 1$: All MCL images have been correctly classified.
- $\text{FPR}_{\text{CLL}} = 0$ and $\text{precision}_{\text{CLL}} = 1$: All CLL images predicted as CLL are correctly classified
- F-measure = 0.946



Class	TPR	TNR	FPR	Precision
CLL	0.857	1	0	1
FL	0.963	0.957	0.043	0.929
MCL	1	0.957	0.043	0.931

Results (3)



Validity of heatmaps

Validity of heatmaps seems to be higher for the proposed ResNet.


Model	Mean [%]	Images without intersection
ResNet50	40	10
Proposed ResNet	60	0

Visual attention mechanism

Accuracy decrease by adding to the model a local branch which focused on a restricted area of the images.



Model	Branch	Training accuracy [%]	Validation accuracy [%]
ResNet50	Global	99.72%	98.67%
	Local	100%	98.67%
	Fusion	100%	98.67%
Proposed ResNet	Global	98.00%	97.33%
	Local	99.00%	83.00%
	Fusion	99.83%	93.33%

Results (4)




Visual attention mechanism

Feeding the local branch with masked images obtained from both good (top) and bad (bottom) intersections → Accuracy decrease




Concluding remarks



- **Artifacts** management improve classification performances.
- High degree of **staining variation** is not a limit for automated lymphoma subtypes classification.
- Splitting images into **patches** is the winning strategy.
- Decreasing **model complexity** leads to slightly lower but comparable classification performances.
- Inclusion of **visual attention mechanism** in the processing pipeline worsen classification performances.
- **Explainable deep learning** could be the meeting point between machine learning and physicians.

Limitations of the study

- Based on a limited dataset including only **3 lymphoma subtypes** and **histological data** (e.g., no molecular data).
- **No clinical confirmation** of the areas of interest indicated by the heatmaps.
 - **Further step:** Ask human hematopathologist to annotate regions of interest in images.



Thanks for your attention!